

Precision and functionality for safe local anaesthesia

PRESTOJECT® Dental Syringes



PRESTOJECT® CW
5374



PRESTOJECT® HK
5376



PRESTOJECT® FL
5375



PRESTOJECT® IL
5373



PRESTOJECT® Safety Stand
5372



Ctra. Sabadell a Granollers, Km. 14,5
Lliçà de Vall, 08185. Barcelona (Spain)
Email: export@inibsa.com
Tel.: +34 938 609 500
Fax: +34 938 437 389

 **CARL MARTIN**
Neuenkamper Straße 80-86
42657 Solingen, Germany
<https://carlmartin.de>

Scan me



60002423 - 01 / 2024

Precision and functionality for safe local anaesthesia

Trust is essential in the dentist-patient relationship and your success is often achieved by offering your patient a comfortable experience. Using safe and precise instruments can help you handle their fear and discomfort.

Prestoject® Dental Syringes are our solution for a successful dental anaesthesia procedure.


PRESTOJECT® DENTAL SYRINGES FEATURES

- Manufactured with stainless steel for long-lasting performance
- Highly resistant to cleaning agents for enhanced longevity
- Ergonomic design for secure and comfortable grip
- Compatible with 1.8ml dental cartridges and metric thread needles
- Various aspiration options available
- Endorsed by Inibsa's expertise in dental pain control management

PRESTOJECT® Dental Syringes

PRESTOJECT® CW

- Cartridge syringe with **craw tip**
- Thumb ring
- For flat plunger cartridges
- Fast and easy side loading
- Manual aspiration
- Suitable for thermal disinfection and 134°C sterilisation



PRESTOJECT® HK

- Cartridge syringe with three **hooks tip**
- Thumb ring
- For hole plunger cartridges
- Folding mechanism loading
- Manual aspiration
- Suitable for thermal disinfection and 134°C sterilisation



PRESTOJECT® FL

- Cartridge syringe with **flat tip**
- Thumb bar
- Compatible with all cartridges with hole or flat plungers
- Fast and easy side loading
- Safe cannula guide into the ampoule
- Self-aspirating: Aspiration is accomplished by easing the pressure on the piston
- Suitable for thermal disinfection and 134° C sterilisation



PRESTOJECT® IL

- Syringe for **intra**ligamental anaesthesia
- Compatible with all cartridges with hole or flat plungers
- Cylinder secured with an easy to use bayonet catch
- Anti-bursting sleeve
- Treble window for optimal control of the anaesthesia cartridge consumption
- Suitable for thermal disinfection and 134°C sterilisation



PRESTOJECT® Safety Stand

- Reduces the risk of injuries or contamination
- Heavy stainless steel for a safer handling
- For safe needle recapping and holding
- For syringe holding during treatment
- Mechanical locking retainer and releasing push button
- Suitable for thermal disinfection and 134° C sterilisation

